



I180803DISB

PromarkerD ELISA vs MS analysis of IGFBP3 concentrations

Scott Bringans, Kirsten Peters

Aim:

Compare IGFBP3 concentrations from ELISA (using new antibodies) and mass spec testing of the same plasma samples using Bland Altman plot analysis. Can the ELISA biomarker concentrations be adjusted to be in agreement with the MS values.

Data:

This analysis was performed on the ELISA plasma data from 2nd August 2018 and the original MS plasma data.

Summary:

- For IGFBP3 there has been a significant improvement in the dynamic range of the ELISA plasma readings with the use of the new antibodies.
- The mean difference between the two tests is 1.92 ug/mL. The ELISA values are higher than the MS values.
- When all ELISA values are decreased by 1.92 ug/mL there is no significant difference between the two methods.




N	98	
	Minimum	Maximum
Plasma IBP3 MS	0.080841619	4.106694445
Plasma rpt IBP3 ELISA (Plasma IBP3 MS + Plasma rpt IBP3 ELISA) / 2	1.896030715	5.546472393
Correlation - r	0.538	


Fit Differences

Parameter	Estimate	95% CI	SE
Mean difference	1.918541907	1.798780376 to 2.038303439	0.060341680
95% Lower LoA	0.747753194	0.542386409 to 0.953119979	0.103473767
95% Upper LoA	3.089330620	2.883963836 to 3.294697405	0.103473767

SD | 0.597352157

- After 1.92 ug/mL was subtracted from all ELISA concentrations, the Bland Altman plot analysis was repeated:





N | 98

	Minimum	Maximum
IGFBP3 MS	0.0808	4.1067
IGFBP3 EL -1.92	-0.0240	3.6265
(IGFBP3 MS + IGFBP3 EL -1.92) / 2	0.0284	3.2751

Correlation - r | 0.538

Fit Differences

Parameter	Estimate	95% CI	SE
Mean difference	-0.0015	-0.12122 to 0.11830	0.06034
95% Lower LoA	-1.1722	-1.37761 to -0.96688	0.10347
95% Upper LoA	1.1693	0.96396 to 1.37470	0.10347

SD | 0.5974

P=0.98 (Methods are not significantly different)

- There were three patients outside the 95% limits of agreement (9027, 9298, and 9176). The following plots remove those three patients.

